

ADULT CRITICAL CARE UNIT BED SPACE SET-UP REFERENCE

ACCU 15TH FLOOR BED LAYOUT



(1) VENTILATOR WITH:
DRY CIRCUIT, BACTERIAL
FILTER
HME FILTER, CO2 SAMPLING

(2) MINDRAY MONITOR
WITH:
ECG, SPO2, NIBP,
DISPOSABLE ADULT CUFF
INVASIVE X2 AND SPLITTER
CO2 MODULE & WATER
TRAP

(3) IV ROLL STAND WITH:
DOCKING STATION
2 VOLUMETRIC PUMPS
2 SYRINGE PUMPS
1 FEEDING PUMP

(4) HIGH VACUUM SUCTION
WITH:
WALL BRACKET
2L SUCTION CANISTER,
VACUUM TUBING, LINER &
VACUUM TUBING

(5) IV ROLL STAND WITH:
DUAL BLOOD PRESSURE
TRANSDUCER PLATE
500ML PRESSURE BAG

OTHER ITEMS:

6 DUAL OXYGEN FLOWMETER & WATERS CIRCUIT

7 LIMB COMPRESSOR AND HOSES

8 CUFF INFLATOR

9 STETHOSCOPE

10 BEDSIDE CHART TROLLEY AND CONTENTS

11 PATIENT BEDSIDE LOCKER AND CONTENTS

12 EMERGENCY AIRWAY BAG (NOT SHOWN)

13 CLINICAL WASTE BIN

14 SHARPS & DRUGS BINS

15 BEDSIDE TABLE

16 TYMPANIC THERMOMETER (INSIDE CHART
TROLLEY)

17 FULL OXYGEN CYLINDER IN CRADLE (NOT
SHOWN)

ADULT CRITICAL CARE UNIT BED SPACE SET-UP REFERENCE



NON INVASIVE
BLOOD PRESSURE
CUFF (NIBP)



NON INVASIVE
BLOOD PRESSURE
CABLE



INVASIVE PRESSURE
SPLITTER
CABLE



INVASIVE PRESSURE
CABLE (x 2)



ECG
CABLE




SP02
(SATS)
CABLE

MINDRAY PATIENT MONITORING

ADULT CRITICAL CARE UNIT BED SPACE SET-UP REFERENCE



LIMB COMPRESSOR
PUMP (FLOWTRON)
2XTUBING



CO2 MODULE
GAS SAMPLING LINE
CO2 WATER TRAP



HME



BACTERIAL FILTER



VENTILATION
DRY CIRCUIT
Part number:
AMVC1812-025

ADULT CRITICAL CARE UNIT BED SPACE SET-UP REFERENCE



TYMPANIC
THERMOMETER



STETHOSCOPE



DOUBLE
TRANSUCER
PLATE



WATERS CIRCUIT



1X 500ML PRESSURE BAG


ADULT CRITICAL CARE UNIT BED SPACE SET-UP REFERENCE



- 1X NG FEED PUMP & POWER CABLE
- 2X VOLUMETRIC PUMP
- 2X SYRINGE PUMP
- 1X PUMP DOCKING STATION & POWER CABLE
- 1X IV POLE



CUFF INFLATOR



- 2 LITRE SUCTION CANISTER
- 2 LITRE SUCTION LINER
- WALL BRACKET CLAMP
- SUCTION TUBING